

Protocol: DNA Purification from a Buccal Brush Using the Gentra Puregene Buccal Cell Kit

This protocol is for purification of genomic DNA from 1 buccal brush using the Gentra Puregene Buccal Cell Kit.

Things to do before starting

- Preheat water baths to 55°C for use in step 3b and 65°C for use in steps 3a and 17 of the procedure.
- Optional: Preheat water bath to 37°C for use in step 5 of the procedure.

Procedure

1. **To collect buccal cells, scrape the inside of the mouth 10 times with a Buccal Collection Brush.**

For best results, wait at least 1 h after eating or drinking to collect buccal cells.

DNA may be purified immediately or samples may be stored on the collection brush for up to 1 month at room temperature (15–25°C).

2. **Dispense 300 µl Cell Lysis Solution into a 1.5 ml microcentrifuge tube. Remove the collection brush from its handle using sterile scissors or a razor blade, and place the detached head in the tube.**

Samples are stable in Cell Lysis Solution for at least 2 years at room temperature.

If 300 µl Cell Lysis Solution is not sufficient to cover the head, the protocol must be scaled up to use a larger volume. Contact QIAGEN Technical Services for more information (see back cover or visit www.qiagen.com).

3. **Complete cell lysis by following step 3a or 3b below:**
 - 3a. **Incubate at 65°C for at least 15 min (up to 60 min for maximum yield).**
 - 3b. **If maximum yield is required, add 1.5 µl Puregene Proteinase K (cat. no. 158918), mix by inverting 25 times, and incubate at 55°C for at least 1 h (up to overnight for maximum yield).**
4. **Remove the collection brush head from the Cell Lysis Solution, scraping it on the sides of the tube to recover as much liquid as possible.**
5. **Optional: If RNA-free DNA is required, add 1.5 µl RNase A Solution, and mix by inverting 25 times. Incubate for 15 min at 37°C. Incubate for 1 min on ice to quickly cool the sample.**

Samples can be incubated at 37°C for up to 1 h.

6. **Add 100 µl Protein Precipitation Solution, and vortex vigorously for 20 s at high speed.**
7. **Incubate for 5 min on ice.**

8. Centrifuge for 3 min at 13,000–16,000 x g.

The precipitated proteins should form a tight pellet. If the protein pellet is not tight, incubate on ice for 5 min and repeat the centrifugation.

9. Pipet 300 µl isopropanol and 0.5 µl Glycogen Solution (cat. no. 158930) into a clean 1.5 ml microcentrifuge tube, and add the supernatant from the previous step by pouring carefully.

Be sure the protein pellet is not dislodged during pouring.

10. Mix by inverting gently 50 times.

11. Centrifuge for 5 min at 13,000–16,000 x g.

12. Carefully discard the supernatant, and drain the tube by inverting on a clean piece of absorbent paper, taking care that the pellet remains in the tube.

13. Add 300 µl of 70% ethanol and invert several times to wash the DNA pellet.

14. Centrifuge for 1 min at 13,000–16,000 x g.

15. Carefully discard the supernatant. Drain the tube on a clean piece of absorbent paper, taking care that the pellet remains in the tube. Allow to air dry for 5 min.

The pellet might be loose and easily dislodged. Avoid over-drying the DNA pellet, as the DNA will be difficult to dissolve.

16. Add 100 µl DNA Hydration Solution and vortex for 5 s at medium speed to mix.

17. Incubate at 65°C for 1 h to dissolve the DNA.

18. Incubate at room temperature overnight with gentle shaking. Ensure tube cap is tightly closed to avoid leakage. Samples can then be centrifuged briefly and transferred to a storage tube.